




Product specification

Module linfusion heater



Product information

Product number	97.8710.550
----------------	-------------

Product descriptions

Name	Module linfusion heater
Product number	97.8710.550
Product description	Module linfusion heater, for warming infusion solutions, tubes, instruments and contrast agent in the SAHARA-III

Size

Width of product	4.5 cm
Length of product	6 cm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Height of case	60 mm
Width of case	130 mm
Depth of case	140 mm
Volume of case	0.0011 cbm

Weight of case	0.1 kg
EAN inner box	4038917968379
EAN case	4038917968386

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT