




Product specification

Collecting tray



Product information

Product number	97.8710.501
----------------	-------------

Product descriptions

Name	Collecting tray
Product number	97.8710.501
Product description	Collecting tray, for collecting plasma leaking from defective storage bags

Size

Width of product	26 cm
Length of product	48 cm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / pallet	72
Height of case	140 mm
Width of case	293 mm
Depth of case	513 mm
Volume of case	0.021 cbm

Weight of case	1.795 kg
EAN inner box	4038917975827
EAN case	4038917975834

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT