



Product specification

SAHARA-III 115 V



Product information

Product number	97.8710.502
----------------	-------------

Product descriptions

Name	SAHARA-III 115 V
Product number	97.8710.502
Product description	SAHARA-III 115 V, for thawing and warming blood products in plastic bags before transfusion, max. capacity: 3 storage bags, suitable for plasma, whole blood, erythrocyte concentrates and cryopreserved preparations

Size

Width of product	32 cm
Length of product	49.3 cm
Height of product	32.5 cm

Product characteristics

Power consumption	655 W
Supply voltage	115 V

Purity & certification

Product category	Medical device equipment, CE
CE certified	CE - manufacturer's self-declaration

Batched	no
---------	----

Packaging

Volume of case	0.129 cbm
Weight of case	16.6 kg
EAN inner box	4038917968331
EAN case	4038917968348

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT