



Product specification

Ink ribbon



Product information

Product number	79.8710.576
----------------	-------------

Product descriptions

Name	Ink ribbon
Product number	79.8710.576
Product description	Ink ribbon, Consumables for protocol printer SP542MD, violet

Material & colours

Color of product	violet
------------------	--------

Size

Width of product	7 cm
Length of product	12 cm

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Piece(s) / pallet	32
Height of case	220 mm
Width of case	390 mm
Depth of case	440 mm

Volume of case	0.0378 cbm
Weight of case	4.5 kg
EAN inner box	4038917975902
EAN case	4038917975919

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT