



Product specification

Ink ribbon



Product information

Product number	79.8710.576
----------------	-------------

Product descriptions

Name	Ink ribbon
Product number	79.8710.576
Product description	Ink ribbon, Consumables for protocol printer SP542MD, violet

Material & colours

Color of product	violet
------------------	--------

Size

Width of product	7 cm
Length of product	12 cm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	144
Piece(s) / pallet	4608
Height of inner box	30 mm

Width of inner box	70 mm
Depth of inner box	120 mm
Height of case	230 mm
Width of case	400 mm
Depth of case	480 mm
Volume of case	0.0442 cbm
Weight of case	7.5 kg
EAN inner box	4038917975902
EAN case	4038917975919

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT