




# Product specification

## Paper roll



### Product information

Product number	79.8710.575
----------------	-------------

### Product descriptions

Name	Paper roll
Product number	79.8710.575
Product description	Paper roll, Consumables for protocol printer, white

### Material & colours

Product material	Paper
Color of product	white

### Size

Width of product	7.5 cm
Length of product	8.3 cm

### Purity & certification

Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order quantity	10
Piece(s) / inner box	10
Piece(s) / case	50

Piece(s) / pallet	3000
Height of case	170 mm
Width of case	400 mm
Depth of case	400 mm
Volume of case	0.0272 cbm
Weight of case	14.15 kg
EAN inner box	4038917975865
EAN case	4038917975872

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT