

Product specification

Mobile blood donation seat



Product information

Product number	97.8700.310
----------------	-------------

Product descriptions

Name	Mobile Blood Donation Seat
------	----------------------------

Product number	97.8700.310
----------------	-------------

Product description	Mobile blood donation seat, for positioning the donor during the blood donation, load bearing capacity: 200 kg, can be collapsed quickly and easily, (LxWxH): 12.5 x 56 x 88 cm
---------------------	---

Size

Width of product	56 cm
------------------	-------

Length of product	12.5 cm
-------------------	---------

Height of product	88 cm
-------------------	-------

Material & colours

Color of product	red
------------------	-----

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched	no
---------	----

Packaging

Weight of case	13 kg
----------------	-------

EAN inner box	4028339671322
---------------	---------------

EAN case	4028339671315
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

