



Product specification

Urine container, 2 l, brown, with light protection, graduated



Product information

Product number	77.580
----------------	--------

Product descriptions

Name	Urine container, 2 l, brown, with light protection, graduated
Product number	77.580
Product description	Urine container, Urine collection, max. working volume: 2 l, Ø opening: 70 mm, brown, with light protection, graduated, material: PE, with paper label, german/english, cap: yellow, Label loosely enclosed, 30 piece(s)/case

Material & colours

Product material	Polyethylene (PE)
Color of product	brown
Closure material	Polypropylene (PP)
Color of cap	yellow

Size

Length excluding cap	235 mm
Length including cap	238 mm
Width of product	106 mm

Length of product	106 mm
Height of product	240 mm
Diameter of opening	70 mm
Max. volume	2 l
Sample volume	2000 ml

Product characteristics

Label/Print	With paper label
Color of print/label	yellow/black
Graduation	yes
Light protection	yes
Viewing window	no
Closure type	Screw cap

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes
Safety data sheet	no

Packaging

Minimum order quantity	30
Piece(s) / inner box	30
Piece(s) / case	30
Piece(s) / pallet	360
Height of case	500 mm
Width of case	343 mm
Depth of case	563 mm
Volume of case	0.0966 cbm
Weight of case	3.5 kg
EAN inner box	4038917082655
EAN case	4038917001700

This document was prepared by EDP support
and is valid without signature.



SARSTEDT