



Product specification

Screw cap, HD-PE, natural, 50 piece(s)/bag



Product information

Product number	76.564
----------------	--------

Product descriptions

Name	Screw cap, HD-PE, natural, 50 piece(s)/bag
Product number	76.564
Product description	Screw cap, suitable for 100-ml container 75.563, material: HD-PE, natural, 50 piece(s)/bag

Material & colours

Product material	High Density Polyethylene (HD-PE)
Color of product	natural

Size

Diameter	62 mm
Height of product	14 mm

Product characteristics

Closure type	Screw cap
Example of use	Cap for urine collection tube

Purity & certification

Product category	diagnostic device, CE
CE certified	CE

Batched	yes
---------	-----

Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	21000
Height of case	250 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0367 cbm
Weight of case	3.91 kg
EAN inner box	4038917082518
EAN case	4038917001595

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT