



## Product specification

**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer**



## Product information

|                |            |
|----------------|------------|
| Product number | 75.562.400 |
|----------------|------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent   |
| Product number      | 75.562.400   |
| Product description | Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: yellow, with safety label, with integrated transfer unit, cap assembled, sterile, 5 piece(s)/tube bag |

## Material & colours

|                  |                                   |
|------------------|-----------------------------------|
| Color of product | transparent                       |
| Closure material | High Density Polyethylene (HD-PE) |
| Color of cap     | yellow                            |

## Size

|                      |       |
|----------------------|-------|
| Diameter             | 57 mm |
| Length excluding cap | 72 mm |

|                     |        |
|---------------------|--------|
| Height of product   | 76 mm  |
| Diameter of opening | 62 mm  |
| Max. volume         | 100 ml |
| Sample volume       | 100 ml |

## Product characteristics

|                      |   |
|----------------------|---|
| Design               | with integrated transfer unit                   |
| Label/Print          | With safety label                               |
| Color of print/label | white/red/black                                 |
| Graduation           | yes   |
| Closure type         | Screw cap                                       |
| Example of use       | closed urine transfer into a V-Monovette® Urine |

## Purity & certification

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Purity standard  | sterile                              |
| Sterilisation    | Electron irradiation                 |
| Batched          | yes                                  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 200           |
| Piece(s) / case      | 200           |
| Piece(s) / pallet    | 5600          |
| Height of case       | 255 mm        |
| Width of case        | 400 mm        |
| Depth of case        | 590 mm        |
| Volume of case       | 0.0602 cbm    |
| Weight of case       | 4.31 kg       |
| EAN inner box        | 4038917236805 |
| EAN case             | 4038917236799 |

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT