




## Product specification

**Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, transparent**



## Product information

Product number	75.562.300
----------------	------------

## Product descriptions

Name	Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, transparent
Product number	75.562.300
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: white, cap enclosed, 50 piece(s)/bag

## Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	white

## Size

Diameter	57 mm
Length excluding cap	72 mm
Height of product	76 mm
Diameter of opening	62 mm

Max. volume	100 ml
Sample volume	100 ml

## Product characteristics

Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	7500
Height of case	330 mm
Width of case	323 mm
Depth of case	443 mm
Volume of case	0.0472 cbm
Weight of case	4.59 kg
EAN inner box	4038917113656
EAN case	4038917274739

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT