



## Product specification

**Urine container, 125 ml, (ØxH): 66 x 67 mm, PP, highly transparent**



## Product information

Product number	75.1355
----------------	---------

## Product descriptions

Name	Urine container, 125 ml, (ØxH): 66 x 67 mm, PP, highly transparent
Product number	75.1355
Product description	Urine container, collection and storage of urine, max. working volume: 125 ml, (ØxH): 66 x 67 mm, Ø opening: 66 mm, highly transparent, graduated, material: PP, 50 piece(s)/bag

## Material & colours

Color of product	highly transparent
------------------	--------------------

## Size

Diameter	66 mm
Length excluding cap	67 mm
Height of product	67 mm
Diameter of opening	66 mm
Max. volume	125 ml
Sample volume	125 ml

## Product characteristics

Graduation	yes
------------	-----

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	24000
Height of inner box	65 mm
Width of inner box	65 mm
Depth of inner box	500 mm
Height of case	265 mm
Width of case	390 mm
Depth of case	590 mm
Volume of case	0.061 cbm
Weight of case	6.2 kg
EAN inner box	4038917409568
EAN case	4042836010088

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT