



Product specification

Urine container, 125 ml, (ØxH): 66 x 67 mm, PP, highly transparent



Product information

Product number	75.1355
----------------	---------

Product descriptions

Name	Urine container, 125 ml, (ØxH): 66 x 67 mm, PP, highly transparent
Product number	75.1355
Product description	Urine container, collection and storage of urine, max. working volume: 125 ml, (ØxH): 66 x 67 mm, Ø opening: 66 mm, highly transparent, graduated, material: PP, 50 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	highly transparent

Size

Diameter	66 mm
Length excluding cap	67 mm
Height of product	67 mm
Diameter of opening	66 mm
Max. volume	125 ml
Sample volume	125 ml

Product characteristics

Graduation	yes
------------	-----

Purity & certification

Product category	diagnostic device, CE
CE certified	CE supplier
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	24000
Height of case	590 mm
Width of case	265 mm
Depth of case	390 mm
Volume of case	0.061 cbm
Weight of case	6.2 kg
EAN inner box	4038917409568
EAN case	4042836010088

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT