



## Product specification

**Multi-purpose container, 125 ml, (LxØ): 85 x 62 mm,  
graduated, PP, transparent**



## Product information

Product number	75.565
----------------	--------

## Product descriptions

Name	Multi-purpose container, 125 ml, (LxØ): 85 x 62 mm, graduated, PP, transparent
Product number	75.565
Product description	Multi-purpose container, max. working volume: 125 ml, (LxØ): 85 x 62 mm, Ø opening: 62 mm, transparent, graduated, material: PP, 50 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Size

Diameter	62 mm
Length excluding cap	85 mm
Length of product	85 mm
Height of product	85 mm
Diameter of opening	62 mm
Max. volume	125 ml

Sample volume	125 ml
---------------	--------

## Product characteristics

Graduation	yes
------------	-----

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	18000
Height of case	195 mm
Width of case	368 mm
Depth of case	578 mm
Volume of case	0.0415 cbm
Weight of case	5.93 kg
EAN inner box	4038917082525
EAN case	4038917001601

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT