



## Product specification

**Container with screw cap, 100 ml, Ø: 57 mm, PP,  
transparent**



## Product information

Product number	75.563
----------------	--------

## Product descriptions

Name	Container with screw cap, 100 ml, Ø: 57 mm, PP, transparent
Product number	75.563
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø: 57 mm, Ø opening: 57 mm, transparent, graduated, material: PP, screw cap, without cap, 50 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Size

Diameter	57 mm
Length excluding cap	72 mm
Length of product	72 mm
Diameter of opening	57 mm
Max. volume	100 ml
Sample volume	100 ml

## Product characteristics

Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	15000
Height of case	330 mm
Width of case	323 mm
Depth of case	443 mm
Volume of case	0.0472 cbm
Weight of case	5.01 kg
EAN inner box	4038917082501
EAN case	4038917001588

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT