



Product specification

Reaction tube, 1.5 ml, PP



Product information

Product number	72.690.999
----------------	------------

Product descriptions

Name	Reaction tube, 1.5 ml, PP
Product number	72.690.999
Product description	Reaction tube, working volume: 1.5 ml, material: PP, colour mix, cap: colour mix, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	Colour mix
Color of cap	Colour mix

Size

Length of product	39 mm
Sample volume	1.5 ml

Product characteristics

Graduation	yes
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	20000 x g

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	5000
Piece(s) / pallet	140000
Height of case	345 mm
Width of case	461 mm
Depth of case	410 mm
Volume of case	0.0649 cbm
Weight of case	5.19 kg
EAN inner box	4038917116756
EAN case	4038917015288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT