



## Product specification

### Reaction tube, 1.5 ml, PP



## Product information

Product number	72.690.999
----------------	------------

## Product descriptions

Name	Reaction tube, 1.5 ml, PP
Product number	72.690.999
Product description	Reaction tube, working volume: 1.5 ml, material: PP, colour mix, cap: colour mix, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

## Material & colours

Color of product	Colour mix
Color of cap	Colour mix

## Size

Length of product	39 mm
Sample volume	1.5 ml

## Product characteristics

Graduation	yes
Cap	Retaining loop cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	20000 x g

## Purity & certification

Batched	yes
---------	-----

## Packaging

Piece(s) / inner box	1200
Piece(s) / case	6000
Height of case	377 mm
Width of case	177 mm
Depth of case	498 mm
Volume of case	0.0332 cbm
Weight of case	8.21 kg
EAN inner box	4038917964920
EAN case	4038917964937

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT