



Product specification

Reaction tube, 1.5 ml, PP



Product information

Product number	72.690.304
----------------	------------

Product descriptions

Name	Reaction tube, 1.5 ml, PP
Product number	72.690.304
Product description	Reaction tube, working volume: 1.5 ml, material: PP, brown, cap: brown, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	brown
Color of cap	brown

Size

Length of product	39 mm
Sample volume	1.5 ml

Product characteristics

Graduation	yes
Cap	Retaining loop cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	20000 x g

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	160000
Volume of case	0.0558 cbm
Weight of case	5.6 kg
EAN inner box	4038917369299
EAN case	4038917369282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT