



Product specification

Reaction tube, 1.5 ml, PP



Product information

Product number	72.690.301
----------------	------------

Product descriptions

Name	Reaction tube, 1.5 ml, PP
Product number	72.690.301
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Color of cap	transparent

Size

Length of product	39 mm
Sample volume	1.5 ml

Product characteristics

Graduation	yes
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	20000 x g

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	160000
Volume of case	0.0558 cbm
Weight of case	5.96 kg
EAN inner box	4038917369268
EAN case	4038917369251

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT