



Product specification

SafeSeal reaction tube, 5 ml, PP, PCR Performance Tested



Product information

Product number	72.701.400
----------------	------------

Product descriptions

Name	SafeSeal reaction tube, 5 ml, PP, PCR Performance Tested
Product number	72.701.400
Product description	SafeSeal reaction tube, working volume: 5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, 100 piece(s)/minigrip bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	17 mm
Length of product	60 mm

Product characteristics

Graduation	yes
------------	-----

Cap	SafeSeal cap
Centrifugation max (RCF)	25000 x g

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	200
Piece(s) / inner box	200
Piece(s) / case	800
Piece(s) / pallet	38400
Height of inner box	124 mm
Width of inner box	184 mm
Depth of inner box	294 mm
Height of case	145 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0345 cbm
Weight of case	3.4628 kg
EAN inner box	4038917398558
EAN case	4038917398541

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT