



## Product specification

### SafeSeal reaction tube, 1.5 ml, PP, PCR Performance Tested, Low DNA-binding



## Product information

Product number	72.706.700
----------------	------------

## Product descriptions

Name	SafeSeal reaction tube, 1.5 ml, PP, PCR Performance Tested, Low DNA-binding
Product number	72.706.700
Product description	SafeSeal reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low DNA-binding, 50 piece(s)/minigrip bag

## Material & colours

Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	transparent

## Size

Diameter	10.8 mm
Length of product	39 mm

## Product characteristics

Graduation	yes
Cap	SafeSeal cap
Centrifugation max (RCF)	30000 x g

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	200
Piece(s) / case	800
Piece(s) / pallet	62400
Height of inner box	96 mm
Width of inner box	150 mm
Depth of inner box	188 mm
Height of case	135 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0188 cbm
Weight of case	1.5104 kg
EAN inner box	4038917350068
EAN case	4038917350051

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT