



Product specification

Reaction tube, 0.5 ml, PP



Product information

Product number	72.698
----------------	--------

Product descriptions

Name	Reaction tube, 0.5 ml, PP
------	---------------------------

Product number	72.698
----------------	--------

Product description	Reaction tube, working volume: 0.5 ml, material: PP, transparent, cap: transparent, cap with ventilation, cap attached, 1,000 piece(s)/bag
---------------------	--

Material & colours

Color of product	transparent
------------------	-------------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	transparent
--------------	-------------

Size

Diameter	7.8 mm
----------	--------

Length of product	30 mm
-------------------	-------

Product characteristics

Graduation	no
------------	----

Cap	cap with ventilation
-----	----------------------

Closure type	Hinged lid
--------------	------------

Centrifugation max (RCF)	20000 x g
--------------------------	-----------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	320000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	2.74 kg
EAN inner box	4038917083508
EAN case	4038917002424

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT