



Product specification

Screw cap micro tube, 0.5 ml



Product information

Product number	72.733.201
----------------	------------

Product descriptions

Name	Screw cap micro tube, 0.5 ml
Product number	72.733.201
Product description	Screw cap micro tube, for volume limitation to 100 µl, working volume: 0.5 ml, conical base, with knurling, transparent, without cap, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	8 mm
Length of product	44 mm

Product characteristics

Base shape	Conical base
Graduation	no
Knurls	yes
Example of use	for volume limitation to 100 µl
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	210000
Height of case	250 mm
Width of case	378 mm
Depth of case	375 mm
Volume of case	0.0354 cbm
Weight of case	4.48 kg
EAN inner box	4038917366625
EAN case	4038917366618

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT