



Product specification

Screw cap micro tube, 1.5 ml, sterile



Product information

Product number	72.687.772
----------------	------------

Product descriptions

Name	Screw cap micro tube, 1.5 ml, sterile
Product number	72.687.772
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, transparent, without cap, sterile, 500 piece(s)/double bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	Conical base
Graduation	no
Knurls	no
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	2500
Piece(s) / inner box	500
Piece(s) / case	2500
Piece(s) / pallet	130000
Height of case	130 mm
Width of case	368 mm
Depth of case	548 mm
Volume of case	0.0262 cbm
Weight of case	3.27 kg
EAN inner box	4038917116329
EAN case	4038917043823

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT