



## Product specification

### Screw cap micro tube, 1.5 ml, PCR Performance Tested



## Product information

Product number	72.692.415
----------------	------------

## Product descriptions

Name	Screw cap micro tube, 1.5 ml, PCR Performance Tested
Product number	72.692.415
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, cap: natural, cap attached assembled, PCR Performance Tested, 100 piece(s)/bag

## Material & colours

Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	Conical base
Graduation	no

Knurls	yes
Cap	Retaining loop cap
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	99000
Height of inner box	141 mm
Width of inner box	191 mm
Depth of inner box	211 mm
Height of case	160 mm
Width of case	223 mm
Depth of case	398 mm
Volume of case	0.0142 cbm
Weight of case	2.09 kg
EAN inner box	4038917240598
EAN case	4038917240581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT