



## Product specification

### Felt marker, black, alcohol-resistant



## Product information

|                |        |
|----------------|--------|
| Product number | 95.950 |
|----------------|--------|

## Product descriptions

|      |                                       |
|------|---------------------------------------|
| Name | Felt marker, black, Alcohol-resistant |
|------|---------------------------------------|

|                |        |
|----------------|--------|
| Product number | 95.950 |
|----------------|--------|

|                     |   |
|---------------------|---|
| Product description | Felt marker, black, alcohol-resistant, application e.g.: for labelling plastic tubes and reaction tubes, 10 piece(s)/case |
|---------------------|---|

## Material & colours

|                  |             |
|------------------|-------------|
| Color of product | white/black |
|------------------|-------------|

## Size

|          |       |
|----------|-------|
| Diameter | 13 mm |
|----------|-------|

|                   |        |
|-------------------|--------|
| Length of product | 133 mm |
|-------------------|--------|

## Product characteristics

|             |            |
|-------------|------------|
| Label/Print | With print |
|-------------|------------|

|                      |       |
|----------------------|-------|
| Color of print/label | black |
|----------------------|-------|

|                |  |
|----------------|--|
| Example of use | for labelling plastic tubes and reaction tubes |
|----------------|--|

## Purity & certification

|                  |               |
|------------------|---------------|
| Product category | not regulated |
|------------------|---------------|

|              |            |
|--------------|------------|
| CE certified | no CE mark |
| Batched      | yes        |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 10            |
| Piece(s) / case      | 500           |
| Piece(s) / pallet    | 45000         |
| Height of inner box  | 27 mm         |
| Width of inner box   | 72 mm         |
| Depth of inner box   | 140 mm        |
| Height of case       | 58 mm         |
| Width of case        | 112 mm        |
| Depth of case        | 145 mm        |
| Volume of case       | 0.0009 cbm    |
| Weight of case       | 9.1 kg        |
| EAN inner box        | 4038917380355 |
| EAN case             | 4038917380348 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT