



## Product specification

**Drop collector, length: 112 mm, PP, tube 75 x 13 mm with funnel cup**



## Product information

Product number	86.290
----------------	--------

## Product descriptions

Name	Drop collector, length: 112 mm, PP, tube 75 x 13 mm with funnel cup
Product number	86.290
Product description	Drop collector, length: 112 mm, material: PP, tube 75 x 13 mm with funnel cup, funnel: HD-PE, for extracting drops for salmonella diagnosis, 250 piece(s)/bag

## Size

Length of product	112 mm
Diameter of opening	40 mm

## Product characteristics

Design	tube 75 x 13 mm with funnel cup
--------	---------------------------------

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Latex-free	yes

Batched	no
Safety data sheet	no

## Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	5.45 kg
EAN inner box	4038917081085
EAN case	4038917035620

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT