



## Product specification

**Stirring rod, length: 120 mm, PP, white**



### Product information

Product number	81.970
----------------	--------

### Product descriptions

Name	Stirring rod, Length: 120 mm, PP, white
Product number	81.970
Product description	Stirring rod, length: 120 mm, material: PP, white, 1,000 piece(s)/bag

### Material & colours

Product material	Polypropylene (PP)
Color of product	white

### Size

Diameter	3 mm
Length of product	120 mm

### Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Batched	no
Safety data sheet	no

### Packaging

Minimum order quantity	10000
Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	400000
Height of case	155 mm
Width of case	268 mm
Depth of case	438 mm
Volume of case	0.0182 cbm
Weight of case	6.23 kg
EAN inner box	4038917269186
EAN case	4038917003506

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT