




Product specification

Filter tip, 300 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill



Product information

Product number	70.3040.375
----------------	-------------

Product descriptions

Name	Filter tip, 300 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill
Product number	70.3040.375
Product description	Filter tip, working volume: 300 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/SingleRefill

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	7 mm
Length of product	60 mm
Max. volume	300 µl

Filter	yes
Graduation	yes

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Height of inner box	148 mm
Width of inner box	105 mm
Depth of inner box	333 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	348 mm
Volume of case	0.0277 cbm
Weight of case	2.8109 kg
EAN inner box	4038917433167
EAN case	4038917433150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT