




Product specification

**Pipette tip, 20 µl, transparent, PCR Performance Tested,
1,000 piece(s)/bag**



Product information

Product number	70.3021
----------------	---------

Product descriptions

Name	Pipette tip, 20 µl, transparent, PCR Performance Tested, 1,000 piece(s)/bag
Product number	70.3021
Product description	Pipette tip, working volume: 20 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit, Abimed, Brand and products of identical design, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Length of product	46 mm
Max. volume	20 µl

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	10000
Piece(s) / inner box	1000
Piece(s) / case	10000
Piece(s) / pallet	720000
Height of case	140 mm
Width of case	263 mm
Depth of case	533 mm
Volume of case	0.0196 cbm
Weight of case	2.12 kg
EAN inner box	4038917525909
EAN case	4038917525893

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT