




Product specification

Gelloader pipette tip, 200 µl, transparent, 96 piece(s)/box



Product information

Product number	70.1190.100
----------------	-------------

Product descriptions

Name	Gelloader pipette tip, 200 µl, transparent, 96 piece(s)/box
Product number	70.1190.100
Product description	Gelloader pipette tip, working volume: 200 µl, transparent, suitable for common pipettors, 96 piece(s)/box

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	67 mm
Max. volume	200 µl

Product characteristics

Filter	no
Graduation	no
Example of use	Gelloader
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	96
Piece(s) / case	960
Piece(s) / pallet	103680
Height of inner box	93 mm
Width of inner box	93 mm
Depth of inner box	125 mm
Height of case	100 mm
Width of case	255 mm
Depth of case	475 mm
Volume of case	0.0121 cbm
Weight of case	1.7299 kg
EAN inner box	4038917344593
EAN case	4038917344586

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT