


Product specification

Gelloader pipette tip, 200 µl, transparent, 96 piece(s)/box



Product information

Order number	70.1190.100
--------------	-------------

Product descriptions

Name	Gelloader pipette tip, 200 µl, transparent, 96 piece(s)/box
Product number	70.1190.100
Product description	Gelloader pipette tip, working volume: 200 µl, transparent, suitable for common pipettors, 96 piece(s)/box

Product characteristics

Filter	no
Graduation	no
Example of use	Gelloader
Autoclavable (max. °C)	121 °C
Suitable for	common pipettors

Size

Length of product	67 mm
Max. volume	200 µl

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Purity & certification

Accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Piece(s) / inner box	96
Piece(s) / case	960
Piece(s) / pallet	103680
Volume of case	0.0126 cbm
Weight of case	1.7299 kg
EAN inner box	4038917344593
EAN case	4038917344586

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

