




Product specification

Filter tip, 200 µl, transparent, Biosphere®, 96 piece(s)/box



Product information

Product number	70.1189.215
----------------	-------------

Product descriptions

Name	Filter tip, 200 µl, transparent, Biosphere®, 96 piece(s)/Box
Product number	70.1189.215
Product description	Filter tip, working volume: 200 µl, transparent, filling level rings, Biosphere®, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Finnpiquette and products of identical design, 96 piece(s)/box

Material & colours

Color of product	transparent
Filter material	Polyethylene (PE)

Size

Length of product	77 mm
Max. volume	200 µl

Product characteristics

Filter	yes
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	Biosphere
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	96
Piece(s) / case	960
Piece(s) / pallet	120960
Height of inner box	85 mm
Width of inner box	84 mm
Depth of inner box	120 mm
Height of case	105 mm
Width of case	250 mm
Depth of case	420 mm
Volume of case	0.011 cbm
Weight of case	1.9402 kg
EAN inner box	4038917314589
EAN case	4038917314572

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT