




Product specification

Pipette tip, 10 ml, transparent, 25 piece(s)/box



Product information

Product number	70.1187.101
----------------	-------------

Product descriptions

Name	Pipette tip, 10 ml, transparent, 25 piece(s)/Box
Product number	70.1187.101
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 25 piece(s)/box

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Length of product	160.8 mm
Max. volume	10 ml

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	250
Piece(s) / inner box	25
Piece(s) / case	250
Piece(s) / pallet	7500
Height of case	610 mm
Width of case	178 mm
Depth of case	305 mm
Volume of case	0.033 cbm
Weight of case	3.63 kg
EAN inner box	4038917126526
EAN case	4038917080422

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT