


## Product specification

**Pipette tip, 5 ml, transparent, 50 piece(s)/box**



### Product information

Product number	70.1183.101
----------------	-------------

### Product descriptions

Name	Pipette tip, 5 ml, transparent, 50 piece(s)/Box
------	---

Product number	70.1183.101
----------------	-------------

Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 50 piece(s)/box
---------------------	---

### Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	transparent
------------------	-------------

### Size

Length of product	123 mm
-------------------	--------

Max. volume	5 ml
-------------	------

### Product characteristics

Filter	no
--------	----

Graduation	yes
------------	-----

Autoclavable (max. °C)	121 °C
------------------------	--------

### Purity & certification

Product category	no medical device   IVD
------------------	-------------------------

Batched

yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	16500
Height of case	610 mm
Width of case	311 mm
Depth of case	165 mm
Volume of case	0.0313 cbm
Weight of case	3.7 kg
EAN inner box	4038917126465
EAN case	4038917080460

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT