




## Product specification

**Pipette tip, 5 ml, transparent, 50 piece(s)/box**



### Product information

Product number	70.1183.101
----------------	-------------

### Product descriptions

Name	Pipette tip, 5 ml, transparent, 50 piece(s)/Box
Product number	70.1183.101
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Gilson, 50 piece(s)/box

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Length of product	123 mm
Max. volume	5 ml

### Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

### Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	yes
---------	-----

## Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	16500
Height of inner box	150 mm
Width of inner box	118 mm
Depth of inner box	147 mm
Height of case	170 mm
Width of case	310 mm
Depth of case	610 mm
Volume of case	0.0321 cbm
Weight of case	3.7 kg
EAN inner box	4038917126465
EAN case	4038917080460

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT