




Product specification

**Filter tip, 1,250 µl, transparent, PCR Performance Tested,
384 piece(s)/bag**



Product information

Product number	70.1186.410
----------------	-------------

Product descriptions

Name	Filter tip, 1,250 µl, transparent, PCR Performance Tested, 384 piece(s)/bag
Product number	70.1186.410
Product description	Filter tip, working volume: 1,250 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 384 piece(s)/bag

Material & colours

Color of product	transparent
Filter material	Polyethylene (PE)

Size

Length of product	103 mm
Max. volume	1250 µl

Product characteristics

Filter	yes
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	384
Piece(s) / case	3840
Piece(s) / pallet	134400
Height of inner box	50 mm
Width of inner box	290 mm
Depth of inner box	400 mm
Height of case	248 mm
Width of case	343 mm
Depth of case	438 mm
Volume of case	0.0373 cbm
Weight of case	4.8192 kg
EAN inner box	4038917126502
EAN case	4038917080361

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT