




Product specification

Pipette tip, 1,000 µl, transparent, 250 piece(s)/bag



Product information

Product number	70.1181
----------------	---------

Product descriptions

Name	Pipette tip, 1,000 µl, transparent, 250 piece(s)/bag
Product number	70.1181
Product description	Pipette tip, working volume: 1,000 µl, transparent, suitable for MLA pipettes, 250 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	86 mm
Max. volume	1000 µl

Product characteristics

Filter	no
Graduation	no
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark

Batched	no
---------	----

Packaging

Piece(s) / inner box	250
Piece(s) / case	5000
Piece(s) / pallet	180000
Height of case	310 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0491 cbm
Weight of case	4.4102 kg
EAN inner box	4038917085656
EAN case	4038917004299

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT