




Product specification

PCR work tray, PC, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes



Product information

Product number	95.987.002
----------------	------------

Product descriptions

Name	PCR work tray, PC, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes
Product number	95.987.002
Product description	PCR work tray, material: PC, transparent, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, 5 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	79 mm
Length of product	115 mm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	100
Piece(s) / pallet	32400
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	1.25 kg
EAN inner box	4038917124157
EAN case	4038917019118

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT