



## Product specification

**Round cuvette, 3.5 ml, height: 51 mm, PP**



### Product information

Product number	68.752
----------------	--------

### Product descriptions

Name	Round cuvette, 3.5 ml, Height: 51 mm, PP
Product number	68.752
Product description	Round cuvette, max. working volume: 3.5 ml, material: PP, height: 51 mm, 1,000 piece(s)/bag

### Size

Diameter	12 mm
Height of product	51 mm
Max. volume	4 ml

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Type of inner packaging	Bag
Piece(s) / inner box	1000
Piece(s) / case	5000

Piece(s) / pallet	100000
Height of case	345 mm
Width of case	405 mm
Depth of case	455 mm
Volume of case	0.0636 cbm
Weight of case	6.61 kg
EAN inner box	4038917083959
EAN case	4038917002844

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT