



Product specification

Round cuvette, 3.5 ml, height: 51 mm, PP



Product information

Product number	68.752
----------------	--------

Product descriptions

Name	Round cuvette, 3.5 ml, Height: 51 mm, PP
Product number	68.752
Product description	Round cuvette, max. working volume: 3.5 ml, material: PP, height: 51 mm, 1,000 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Size

Diameter	12 mm
Height of product	51 mm
Max. volume	4 ml

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Minimum order quantity	5000
Type of inner packaging	Bag
Piece(s) / inner box	1000
Piece(s) / case	5000
Piece(s) / pallet	100000
Height of case	350 mm
Width of case	408 mm
Depth of case	458 mm
Volume of case	0.0654 cbm
Weight of case	6.61 kg
EAN inner box	4038917083959
EAN case	4038917002844

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT