



Product specification

**Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm, PS,
transparent, optical sides: 2**



Product information

Product number	67.746
----------------	--------

Product descriptions

Name	Semi-micro cuvette, 3 ml, (HxW): 45 x 12 mm, PS, transparent, optical sides: 2
Product number	67.746
Product description	Semi-micro cuvette, max. working volume: 3 ml, material: PS, transparent, (HxW): 45 x 12 mm, optical sides: 2, path length: 10 mm, 500 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Width of product	12 mm
Height of product	45 mm
Max. volume	3 ml

Product characteristics

Path length	10 mm
-------------	-------

Wavelength in the cuvette area

330 nm

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	2000
Type of inner packaging	Bag
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	96000
Height of case	215 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0315 cbm
Weight of case	5.02 kg
EAN inner box	4038917083928
EAN case	4038917002806

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT