




## Product specification

**PCR plate without skirt, 96 well, white, Low Profile, 100 µl,  
PCR Performance Tested, PP**



## Product information

Product number	72.1977.232
----------------	-------------

## Product descriptions

Name	PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP
Product number	72.1977.232
Product description	PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 20 piece(s)/minigrip bag

## Material & colours

Color of product	white
Profile	Low Profile

## Size

Width of product	77.7 mm
Length of product	117.7 mm
Max. volume	100 µl
Sample volume	0.1 ml

Number of wells	96
Construction of the edges of the plate	without skirt
Base shape	conical
Closure type	Film
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	20
Piece(s) / case	100
Piece(s) / pallet	5000
Height of inner box	100 mm
Width of inner box	140 mm
Depth of inner box	155 mm
Height of case	157 mm
Width of case	163 mm
Depth of case	510 mm
Volume of case	0.0131 cbm
Weight of case	1.5 kg
EAN inner box	4038917405133
EAN case	4038917405126

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

