




Product specification

**PCR plate without skirt, 96 well, white, Low Profile, 100 µl,
PCR Performance Tested, PP**



Product information

Product number	72.1977.232
----------------	-------------

Product descriptions

Name	PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP
Product number	72.1977.232
Product description	PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 20 piece(s)/minigrip bag

Material & colours

Product material	Polypropylene (PP)
Color of product	white
Profile	Low Profile

Size

Width of product	77.7 mm
Length of product	117.7 mm
Max. volume	100 µl
Sample volume	0.1 ml

Product characteristics

Number of wells	96
Construction of the edges of the plate	without skirt
Base shape	conical
Closure type	Film
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	100
Piece(s) / inner box	20
Piece(s) / case	100
Piece(s) / pallet	5000
Height of case	493 mm
Width of case	137 mm
Depth of case	137 mm
Volume of case	0.0093 cbm
Weight of case	1.5 kg
EAN inner box	4038917405133
EAN case	4038917405126

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT

This document was prepared by EDP support and is valid without signature.