



Product specification

**PCR plate full skirt, 384 well, white, Low Profile, 40 µl,
PCR Performance Tested, PP**



Product information

| | |
|----------------|-------------|
| Product number | 72.1985.202 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|---|
| Name | PCR plate full skirt, 384 well, white, Low Profile, 40 µl, PCR Performance Tested, PP |
| Product number | 72.1985.202 |
| Product description | PCR plate full skirt, for LightCycler® 480, 384 well, white, Low Profile, 40 µl, PCR Performance Tested, material: PP, optimised for qPCR, 50 piece(s)/minigrip bag |

Material & colours

| | |
|------------------|-------------|
| Color of product | white |
| Profile | Low Profile |

Size

| | |
|-------------------|----------|
| Diameter | 3.3 mm |
| Width of product | 85.5 mm |
| Length of product | 127.8 mm |
| Height of product | 9.6 mm |
| Max. volume | 40 µl |
| Sample volume | 40 µl |

Product characteristics

| | |
|--|------------|
| Number of wells | 384 |
| Construction of the edges of the plate | full skirt |
| Base shape | conical |
| Autoclavable (max. °C) | 121 °C |
| Suitable for qPCR | yes |
| Centrifugation max (RCF) | 3000 x g |
| Centrifugation temperature | 22 °C |
| Centrifugation time | 5 Min |

Purity & certification

| | |
|------------------|------------------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 50 |
| Piece(s) / case | 100 |
| Piece(s) / pallet | 7200 |
| Height of inner box | 138 mm |
| Width of inner box | 180 mm |
| Depth of inner box | 210 mm |
| Height of case | 191 mm |
| Width of case | 270 mm |
| Depth of case | 327 mm |
| Volume of case | 0.0169 cbm |
| Weight of case | 2.15 kg |
| EAN inner box | 4038917251259 |
| EAN case | 4038917251242 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT