



Product specification

PCR plate full skirt, 384 well, transparent, Low Profile, 40 µl, PCR Performance Tested, PP



Product information

Product number	72.1984.202
----------------	-------------

Product descriptions

Name	PCR plate full skirt, 384 well, transparent, Low Profile, 40 µl, PCR Performance Tested, PP
Product number	72.1984.202
Product description	PCR plate full skirt, 384 well, transparent, Low Profile, 40 µl, PCR Performance Tested, material: PP, 50 piece(s)/minigrip bag

Material & colours

Color of product	transparent
Profile	Low Profile

Size

Diameter	3.2 mm
Width of product	85.5 mm
Length of product	127.8 mm
Height of product	9.65 mm
Max. volume	40 µl
Sample volume	40 µl

Product characteristics

Number of wells	384
Construction of the edges of the plate	full skirt
Base shape	conical
Closure type	Film
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	5050
Height of inner box	254 mm
Width of inner box	152 mm
Depth of inner box	203 mm
Height of case	254 mm
Width of case	152 mm
Depth of case	203 mm
Volume of case	0.0078 cbm
Weight of case	2.18 kg
EAN inner box	4038917135016
EAN case	4038917067584

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT