



## Product specification

**PCR plate half skirt, 96 well, white, Low Profile, 100 µl,  
PCR Performance Tested, PP**



## Product information

Product number	72.1982.202
----------------	-------------

## Product descriptions

Name	PCR plate half skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP
Product number	72.1982.202
Product description	PCR plate half skirt, for LightCycler® 480, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 25 piece(s)/minigrip bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	white
Profile	Low Profile

## Size

Diameter	5.38 mm
Width of product	85.3 mm
Length of product	127.8 mm
Height of product	15.55 mm
Max. volume	100 µl

Sample volume

0.1 ml

## Product characteristics

Number of wells	96
Construction of the edges of the plate	half skirt
Base shape	conical
Closure type	Film
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

## Purity & certification

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Minimum order quantity	100
Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	7200
Height of case	335 mm
Width of case	208 mm
Depth of case	278 mm
Volume of case	0.0194 cbm
Weight of case	2.3 kg
EAN inner box	4038917282017
EAN case	4038917282000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

