




Product specification

**PCR strip of lids, PCR Performance Tested, PP,
transparent, flat cap**



Product information

Product number	65.1998.400
----------------	-------------

Product descriptions

Name	PCR strip of lids, PCR Performance Tested, PP, transparent, flat cap
Product number	65.1998.400
Product description	PCR strip of lids, for PCR plates (REF 72.1979.102/.132, 72.1982.202), PCR Performance Tested, material: PP, transparent, flat cap, suitable for qPCR, 12 piece(s)/minigrip bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Width of product	7.4 mm
Length of product	80.5 mm

Product characteristics

Cap	flat cap
Closure type	Push cap

Autoclavable (max. °C)	121 °C
Pierceable	no
Suitable for qPCR	yes

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	1200
Piece(s) / inner box	300
Piece(s) / case	1200
Piece(s) / pallet	151200
Height of case	123 mm
Width of case	255 mm
Depth of case	355 mm
Volume of case	0.0111 cbm
Weight of case	0.69 kg
EAN inner box	4038917336369
EAN case	4038917336352

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT