




## Product specification

**PCR strip of 8, 100 µl, PCR Performance Tested,  
transparent, PP, flat cap**



## Product information

Product number	72.991.103
----------------	------------

## Product descriptions

Name	PCR strip of 8, 100 µl, PCR Performance Tested, transparent, PP, flat cap
Product number	72.991.103
Product description	PCR strip of 8, up to 100 µl, Low Profile, PCR Performance Tested, transparent, material: PP, flat cap, cap attached, suitable for qPCR, 12 piece(s)/minigrip bag

## Material & colours

Color of product	transparent
Profile	Low Profile

## Size

Diameter	6.26 mm
Width of product	10.6 mm
Length of product	71.6 mm
Height of product	15.4 mm
Max. volume	100 µl
Sample volume	0.1 ml

## Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	120
Piece(s) / case	1200
Piece(s) / pallet	43200
Height of inner box	100 mm
Width of inner box	140 mm
Depth of inner box	155 mm
Height of case	305 mm
Width of case	279 mm
Depth of case	305 mm
Volume of case	0.026 cbm
Weight of case	2.1996 kg
EAN inner box	4038917231688
EAN case	4038917231671

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT