




# Product specification

**PCR strip of 8, 200 µl, PCR Performance Tested, colour mix, PP**



## Product information

Product number	72.985.992
----------------	------------

## Product descriptions

Name	PCR strip of 8, 200 µl, PCR Performance Tested, colour mix, PP
Product number	72.985.992
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, red/green/blue/violet, material: PP, without cap, suitable for qPCR, per colour, 120 piece(s)/minigrip bag

## Material & colours

Color of product	Colour mix
Profile	High Profile

## Size

Diameter	5.4 mm
Width of product	8.1 mm
Length of product	70.9 mm
Height of product	20.7 mm
Max. volume	200 µl
Sample volume	0.2 ml

## Product characteristics

Base shape	conical
Graduation	no
Closure type	Push cap
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	480
Piece(s) / pallet	155520
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.6398 kg
EAN inner box	4038917124119
EAN case	4038917067683

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT